

# SSX2IP Antibody



PACO12483

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SSX2IP

**Uniprot**

Q9Y2D8

**Synonyms:**

synovial sarcoma, X breakpoint 2 interacting protein; SSX2IP; ADIP; FLJ10848; KIAA0923; MGC75026 ;

**Immunogen:**

Human SSX2IP.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A