

# SSX4 Antibody



PACO12484

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SSX4

**Uniprot**

O60224

**Synonyms:**

synovial sarcoma, X breakpoint 4; SSX4; MGC119056; MGC12411 ;

**Immunogen:**

Human SSX4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A