

ST3GAL4 Antibody



PACO12487

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ST3GAL4

Uniprot

Q11206

Synonyms:

ST3 beta-galactoside alpha-2,3-sialyltransferase 4; ST3GAL4; CGS23; FLJ11867; FLJ46764; NANTA3; SAT3; SIAT4; SIAT4C; ST3GalIV; STZ ;

Immunogen:

Human ST3GAL4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A