

ST3GAL6 Antibody



PACO12489

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ST3GAL6

Uniprot

Q9Y274

Synonyms:

ST3 beta-galactoside alpha-2,3-sialyltransferase 6; ST3GAL6; SIAT10; ST3GALVI ;

Immunogen:

Human ST3GAL6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A