

ST5 Antibody



PACO12490

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ST5

Uniprot

P78524

Synonyms:

suppression of tumorigenicity 5; ST5; DENND2B; HTS1; MGC33090; p126 ;

Immunogen:

Human ST5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A