

ST6GALNAC1 Antibody



PACO12492

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ST6GALNAC1

Uniprot

Q9NSC7

Synonyms:

ST6 (alpha-N-acetyl-neuraminy-2,3-beta-galactosyl-1,3)-N-acetylgalactosaminide alpha-2,6-sialyltransferase 1; ST6GALNAC1; HSY11339; SIAT7A; ST6GalNAc; STYI ;

Immunogen:

Human ST6GALNAC1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A