

ST7 Antibody



PACO12494

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ST7

Uniprot

Q9NRC1

Synonyms:

suppression of tumorigenicity 7; ST7; DKFZp762O2113; ETS7q; FAM4A; FAM4A1; HELG; RAY1; SEN4; TSG7 ;

Immunogen:

Human ST7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A