

STAM Antibody



PACO12501

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STAM

Uniprot

Q92783

Synonyms:

signal transducing adaptor molecule (SH3 domain and ITAM motif) 1; STAM;
DKFZp686J2352; STAM1 ;

Immunogen:

Human STAM.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A