

# STARD10 Antibody



PACO12505

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

STARD10

**Uniprot**

Q9Y365

**Synonyms:**

StAR-related lipid transfer (START) domain containing 10; STARD10; CGI-52; MGC14401; NY-CO-28; PCTP2; SDCCAG28 ;

**Immunogen:**

Human STARD10.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A