

STAT2 Antibody



PACO12511

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

STAT2

Uniprot

P52630

Synonyms:

signal transducer and activator of transcription 2, 113kDa; STAT2; ISGF-3; MGC59816; P113; STAT113 ;

Immunogen:

Human STAT2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A