

STAT4 Antibody



PACO12513

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STAT4

Uniprot

Q14765

Synonyms:

signal transducer and activator of transcription 4; STAT4; SLEB11 ;

Immunogen:

Human STAT4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A