

STAT6 Antibody



PACO12516

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

STAT6

Uniprot

P42226

Synonyms:

signal transducer and activator of transcription 6, interleukin-4 induced; STAT6; D12S1644; IL-4-STAT; STAT6B; STAT6C ;

Immunogen:

Human STAT6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A