

STIM1 Antibody



PACO12522

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

STIM1

Uniprot

Q13586

Synonyms:

stromal interaction molecule 1; STIM1; D11S4896E; GOK ;

Immunogen:

Human STIM1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A