

STK11IP Antibody



PACO12527

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

STK11IP

Uniprot

Q8N1F8

Synonyms:

serine/threonine kinase 11 interacting protein; STK11IP; KIAA1898; LIP1; LKB1IP; STK11IP1 ;

Immunogen:

Human STK11IP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A