

STK3 Antibody



PACO12531

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

STK3

Uniprot

Q13188

Synonyms:

serine/threonine kinase 3 (STE20 homolog, yeast); STK3; FLJ90748; KRS1; MST2 ;

Immunogen:

Human STK3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A