

STK32B Antibody



PACO12533

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STK32B

Uniprot

Q9NY57

Synonyms:

serine/threonine kinase 32B; STK32B; HSA250839; STK32; STKG6; YANK2 ;

Immunogen:

Human STK32B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A