

STK38L Antibody



PACO12538

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

STK38L

Uniprot

Q9Y2H1

Synonyms:

serine/threonine kinase 38 like; STK38L; KIAA0965; NDR2 ;

Immunogen:

Human STK38L.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A