

STOML1 Antibody



PACO12543

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

STOML1

Uniprot

Q9UBI4

Synonyms:

stomatin (EPB72)-like 1; STOML1; FLJ36370; SLP-1; STORP; hUNC-24 ;

Immunogen:

Human STOML1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A