

STRAP Antibody



PACO12547

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

STRAP

Uniprot

Q9Y3F4

Synonyms:

serine/threonine kinase receptor associated protein; STRAP; MAWD; PT-WD; UNRIP ;

Immunogen:

Human STRAP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A