

STT3A Antibody



PACO12551

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

STT3A

Uniprot

P46977

Synonyms:

STT3, subunit of the oligosaccharyltransferase complex, homolog A (S. cerevisiae); STT3A; FLJ27038; ITM1; MGC9042; STT3-A; TMC ;

Immunogen:

Human STT3A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A