

STT3B Antibody



PACO12552

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STT3B

Uniprot

Q8TCJ2

Synonyms:

STT3, subunit of the oligosaccharyltransferase complex, homolog B (S. cerevisiae); STT3B; FLJ90106; SIMP; STT3-B ;

Immunogen:

Human STT3B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A