

STX11 Antibody



PACO12554

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

STX11

Uniprot

O75558

Synonyms:

syntaxin 11; STX11; FHL4; HLH4; HPLH4 ; Syntaxin-11, SYX11, syntaxin11

Immunogen:

Human STX11.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A