

STXBP6 Antibody



PACO12568

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

STXBP6

Uniprot

Q8NFX7

Synonyms:

syntaxin binding protein 6 (amisyn); STXBP6; FLJ39638; HSPC156; amisyn ;

Immunogen:

Human STXBP6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A