

# SUCLG1 Antibody



PACO12572

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SUCLG1

**Uniprot**

P53597

**Synonyms:**

succinate-CoA ligase, alpha subunit; SUCLG1; FLJ21114; G-ALPHA; SUCLA1 ;

**Immunogen:**

Human SUCLG1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A