

# SUPT16H Antibody



PACO12587

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SUPT16H

**Uniprot**

Q9Y5B9

**Synonyms:**

suppressor of Ty 16 homolog (S. cerevisiae); SUPT16H; CDC68; FACTP140; FLJ10857; FLJ14010; FLJ34357; SPT16/CDC68 ;

**Immunogen:**

Human SUPT16H.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A