

SV2C Antibody



PACO12594

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SV2C

Uniprot

Q496J9

Synonyms:

synaptic vesicle glycoprotein 2C; SV2C; MGC118874; MGC118875; MGC118876; MGC118877 ;

Immunogen:

Human SV2C.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A