

TACC1 Antibody



PACO12618

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TACC1

Uniprot

O75410

Synonyms:

transforming, acid, c coiled-coil containing protein 1; TACC1; DKFZp686K18126; FLJ42304; Ga55; KIAA1103 ;

Immunogen:

Human TACC1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A