

TAF1C Antibody



PACO12625

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TAF1C

Uniprot

Q15572

Synonyms:

TATA box binding protein (TBP)-associated factor, RNA polymerase I, C, 110kDa; TAF1C; MGC:39976; SL1; TAFI110; TAFI95 ;

Immunogen:

Human TAF1C.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A