

# TAGLN3 Antibody



PACO12634

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TAGLN3

**Uniprot**

Q9UI15

**Synonyms:**

199aa,22.4kd;

**Immunogen:**

Human TAGLN3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.320&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A