

TARDBP Antibody



PACO12639

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat, Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF, IP

Recommended dilutions:

Protein Background:

Gene ID:

TARDBP

Uniprot

Q13148

Synonyms:

TAR DNA binding protein; TARDBP; ALS10; TDP-43; Tdp43 ; TDP 43

Immunogen:

Human TARDBP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A