

TARS Antibody



PACO12640

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TARS

Uniprot

P26639

Synonyms:

threonyl-tRNA synthetase; TARS; MGC9344; ThrRS ;

Immunogen:

Human TARS.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A