

# TBC1D5 Antibody



PACO12646

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TBC1D5

**Uniprot**

Q92609

**Synonyms:**

TBC1 domain family, member 5; TBC1D5; KIAA0210 ;

**Immunogen:**

Human TBC1D5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A