

TBRG1 Antibody



PACO12657

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TBRG1

Uniprot

Q3YBR2

Synonyms:

transforming growth factor beta regulator 1; TBRG1; FLJ14621; FLJ25020; FLJ90113; MGC129890; NIAM; TB-5 ;

Immunogen:

Human TBRG1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A