

TBX3 Antibody



PACO12661

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TBX3

Uniprot

O15119

Synonyms:

T-box 3; TBX3; TBX3-ISO; UMS; XHL ;

Immunogen:

Human TBX3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A