

TBX6 Antibody



PACO12663

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TBX6

Uniprot

O95947

Synonyms:

T-box 6; TBX6; DFN67 ;

Immunogen:

Human TBX6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A