

TBXAS1 Antibody



PACO12664

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TBXAS1

Uniprot

P24557

Synonyms:

thromboxane A synthase 1 (platelet); TBXAS1; CYP5; CYP5A1; GHOSAL; THAS; TS; TXAS; TXS ;

Immunogen:

Human TBXAS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A