

TCEAL7 Antibody



PACO12667

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TCEAL7

Uniprot

Q9BRU2

Synonyms:

transcription elongation factor A (SII)-like 7; TCEAL7; MGC23947; MPMGp800C04260Q003 ;

Immunogen:

Human TCEAL7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A