

TCEB3 Antibody



PACO12670

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ELOA

Uniprot

Q14241

Synonyms:

transcription elongation factor B (SIII), polypeptide 3 (110kDa, elongin A); TCEB3; FLJ38760; FLJ42849; SIII; TCEB3A ;

Immunogen:

Human TCEB3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A