

TCP11 Antibody



PACO12682

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TCP11

Uniprot

Q8WWU5

Synonyms:

t-complex 11 homolog (mouse); TCP11; D6S230E; KIAA0229; MGC111103 ;

Immunogen:

Human TCP11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A