

TEK Antibody



PACO12693

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TEK

Uniprot

Q02763

Synonyms:

TEK tyrosine kinase, endothelial; TEK; CD202B; TIE-2; TIE2; VMCM; VMCM1 ;

Immunogen:

Human TEK.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A