

TERF2IP Antibody



PACO12697

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TERF2IP

Uniprot

Q9NYB0

Synonyms:

telomeric repeat binding factor 2, interacting protein; TERF2IP; DRIP5; RAP1 ;

Immunogen:

Human TERF2IP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A