

# TFCP2 Antibody



PACO12711

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TFCP2

**Uniprot**

Q12800

**Synonyms:**

transcription factor CP2; TFCP2; CP2; LBP-1C; LSF; SEF; TFCP2C ;

**Immunogen:**

Human TFCP2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A