

TFDP1 Antibody



PACO12712

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TFDP1

Uniprot

Q14186

Synonyms:

transcription factor Dp-1; TFDP1; DP1; DRTF1; Dp-1 ;

Immunogen:

Human TFDP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A