

TFEB Antibody



PACO12715

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IP, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TFEB

Uniprot

P19484

Synonyms:

transcription factor EB; TFEB; AlphaTFEB; TCFEB; bHLHe35 ;

Immunogen:

Human TFEB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A