

TGFB11 Antibody



PACO12725

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TGFB11

Uniprot

O43294

Synonyms:

transforming growth factor beta 1 induced transcript 1; TGFB11; ARA55; HIC-5; HIC5; TSC-5 ;

Immunogen:

Human TGFB11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A