

TGM1 Antibody



PACO12728

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TGM1

Uniprot

P22735

Synonyms:

transglutaminase 1 (K polypeptide epidermal type I, protein-glutamine-gamma-glutamyltransferase); TGM1; ICR2; KTG; LI; LI1; TGASE; TGK ;

Immunogen:

Human TGM1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A