

# TGM2 Antibody



PACO12729

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TGM2

**Uniprot**

P21980

**Synonyms:**

transglutaminase 2 (C polypeptide, protein-glutamine-gamma-glutamyltransferase); TGM2; G-ALPHA-h; GNAH; TG2; TGC ;

**Immunogen:**

Human TGM2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A