

TGS1 Antibody



PACO12731

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TGS1

Uniprot

Q96RS0

Synonyms:

trimethylguanosine synthase homolog (S. cerevisiae); TGS1; DKFZp762A163; FLJ22995; NCOA6IP; PIMT; PIPMT ;

Immunogen:

Human TGS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A