

THADA Antibody



PACO12733

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

THADA

Uniprot

Q6YHU6

Synonyms:

thyroid adenoma associated; THADA; FLJ21877; FLJ44016; FLJ44876; FLJ77530; GITA; KIAA1767 ;

Immunogen:

Human THADA.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A