

# ACOT13 Antibody



PACO12737

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ACOT13

**Uniprot**

Q9NPJ3

**Synonyms:**

thioesterase superfamily member 2; THEM2; HT012; MGC4961; PNAS-27 ;

**Immunogen:**

Human THEM2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A